Vidrio Industrial POBEL, S.A.

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

1.2. General information

- 1) Handle the parcel with care. Unpack and check that the contents coincide with the packing-list. If any part is damaged or missing, please, notify to your supplier immediately.
- 2) Do not install or use the equipment without reading this handbook before.
- 3) This handbook must always be attached to the equipment and it must be available for all users.
- 4) **IMPORTANT! WE WILL NOT ACCEPT ANY APPARATUS TO BE REPAIRED IF IT IS NOT DULY CLEANED**
- 5) If any modification, elimination or lack of maintenance of any device of the equipment by the user transgress the directive 89/655/CEE, the manufacturer is not responsible for the damage that can occur.
- 6) Do not use the apparatus with liquids which can give off vapours capable of making explosive mixtures.

1.3. Packing list

The standard equipment consists of the following components:

Description	Qty.
Distiller	1
Condenser*	1
Heater	1
Filter	1
Feed tube	1
Support for water still**	1
Set of silicone tubes	1

*Filled with salt for protection

**Built-in with the distiller

1.4. Recommendations

- Electric:
 - Electrical wiring must be, at least, of 2.5 mm² section.
 - It must have ground.
 - Power point must let intensity of 13 A pass.
- Outlet:
 - The tap water should have screw mouth of 3/4".

1.5. Technical features

<i>Code</i>	<i>Length (cm)</i>	<i>Width (cm)</i>	<i>Height (cm)</i>	<i>Weight (Kg)</i>
710 710	55	15	50	4

<i>Code</i>	<i>Output (lts/h)</i>	<i>Conductivity at 20 °C</i>	<i>Water consumption (l/min)</i>	<i>Consumption (W)</i>	<i>Voltage (V)</i>	<i>Intensity (A)</i>
710 710	4	2 μ S / cm	1,0	1 x 3000	230	13


Fig. 1


Fig. 2

1.6.Assembly

Steps to follow:

1. Be sure that heater elements (washer of silicone rubber, guide fitting and PBTB nut) are in the order and position that the fig. 1 shows.


2. Insert the heater through the screw mouth (see fig. 2); hold the heater of the fitting (fig. 1, piece 3) and turn the PTPB nut (fig. 1, piece 2)

3. Place the condenser onto the vertical tube. Be sure that the O-rings are on the surface of the tube and connect the silicone tubes.

4. Put the water inlet tube on the inlet still (fig. 2, piece 6)

5. Be sure of putting the filter, of tightening and connecting the other end to the tap.

- Once you have assembled, check:
 - The silicone tubes hold suitable position to let water pass and to avoid water pressures.

1.7. Operating

Steps to follow:

1. Connect the end wiring of the heater to the socket, which is on the metallic base (piece 12).
2. Turn the tap on.
3. Turn the tap off when you see that the water level is on the heater.

1.8. Switch on

Steps to follow:


1. Switch on the green switch (piece 10). This will light showing that the still is connecting to the power point. The warning light will also light (piece 9) showing that the water supply is not enough.

2. Turn the tap on slowly until the pilot (piece 9) turn off.

- Before few minutes, the heater will be getting red and the boiling will start. The stilled water goes out through the outlet (piece 4).

- If water steam went out through steam hole (piece 2), you would have to increase water supply.

- You should clean the still regularly. So, in this way, the results will be better and you can get a still more durable.



1.9.Maintenance

- You should clean the still regularly. So, in this way, the results will be better and you can get a still more durable.
- The cleaning periodicity depends on the salts which water supply contains.
- You should clean the filters of the tubes whenever you clean the still.

Steps to follow:

1. Put acetic acid, hydrochloric acid or whichever another product into the boiler through cleaning mouth (piece 5). The water should be hot.
2. You have to wait for the attack to the residues.
3. Rinse with clean water once or twice.
4. Fill with water until you see the water level is on the heater.

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VIDRIO INDUSTRIAL POBEL, S.A.