**Vidrio Industrial POBEL, S.A.**

**Destilador de agua / *Water Still***

**DESA 0041 / DESA 0081**

**MANUAL DE INSTRUCCIONES  
*USING MANUAL***

C/ Luis I, s/n - Almacén 12 - Pol. Ind. Vallecas IV  
28031 Madrid (ESPAÑA)  
Tel. (34) 91 380 33 18 - Fax. (34) 91 380 32 16  
E-mail: [pobel@pobel.com](mailto:pobel@pobel.com) - Web: [www.pobel.com](http://www.pobel.com)

## 1.1. Tabla de contenidos

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## 1.2. Getting started

*Thank you very much for purchasing POBEL product.*

*Your water still has been designed with function, reliability and safety in mind. It is your responsibility to install it in conformance with local electrical codes. For safe operation, please pay attention to the alert symbols through the manual.*

*This manual contains important operating and safety information. You must carefully read and understand the contents of this manual prior to the use of this equipment.*

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### Warning

Warning alert you to a possibility of personal injury.

### Caution

Caution alert you to a possibility of damage to the equipment.

### Note

Notes alert you to pertinent facts and conditions.

### Hot surface

Hot surface alert you possibility of burning injury by hot instrument surface.

### Explosive

Explosive alerts you to possibility of explosion by high pressure.

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
### 1.3. Product overview

*POBEL water still produces distilled water from tap water with space saving design and effective distilled water production.*

*Unique design and safety feature with hydraulic temperature controller and a water float switch provides automatic operation. By lack of water, the heating is automatically switched off to protect heater from over heating.*

*Water still is easy to clean-up by just removing upper part of the body (condenser unit). By controlling cooling water flow, you can save a lot of water wasted before.*

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## 1.6. Operating

### ● Before operating

1) The main voltage must correspond to the voltage given on the name-plate.

2) Electrical wiring must be, at least,

	<u>110 V</u>	<u>230 V</u>
DESA 1041	22 Amp.	
DESA 0041	-----	14 Amp.
DESA 0081	-----	26 Amp.

3) Place water still on a flat surface.

4) Place water still near by water supply

### ● Tubing connection

1) Connect the water inlet tap (4) with a hose to the water main.

2) Connect a hose to the water outlet (3) and place its end in a drain on a lower level.

3) Connect a hose to the distilled water outlet (13) and place its end in a distilled water reservoir on a lower level.

### ● Operation

1) Turn the water supply on and check if there is any water leakage between the tubing connections.

2) Switch on (15). The pilot lamp in the main switch will glow.

3) The water still will automatically starts when water in the chamber reaches appropriate level.

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## 1.7. Warning


1. *Make sure the electric supply is in accordance to the given voltage on the name-plate, 50 / 60 Hz, single phase.*
2. *Place water still on a flat surface.*
3. *Do not touch top of the instrument. During operation, surface of the top of the water still and steam exhausted from the hole is very hot.*
4. *Lime deposit should be removed from the evaporator once or twice a month depending on the degree of hardness of the tap water. Take off the upper part of the water still and remove the silicone packing. A mixture of 10% formic acid, 10% acetic acid and 80% distilled water is recommended for cleaning the apparatus. Fill the solvent into the evaporator to just above the highest lime deposit and switch on the water still. The evaporator is clean after approx 10 to 15 minutes. Pour out the lime and solvent mixture thoroughly rise with water.*
5. *In the event of overheating, it activates the security device and you should reset the water still pushing RESET (16).*

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<b>Trouble</b>	<b>Check first</b>	<b>Trouble Shooting</b>
Power failure	Check Electric Supply	Plug firmly into the electric supply
Power on but water do not boil	Check water supply	Supply water enough and appropriate pressure
	Check heater	Replace heater
Contact sales representative or customer service department		


## 1.9. Diagrama de flujo de agua / *Water flow diagram*


## 1.10. Diagrama del circuito / *Circuit diagram*

- 1 -> Interruptor / *Switch*
- 2 -> Termostato seguridad / *Security thermostat*
- 3 -> Control nivel deposito / *Level controller*
- C1 -> Contactor
- C2 -> Contactor
- R1 -> Resistencia / *Heater*
- R2 -> Resistencia / *Heater*

**DESA 0081**

- 1 -> Interruptor / *Switch*
- 2 -> Termostato seguridad / *Security thermostat*
- 3 -> Control nivel deposito / *Level controller*
- C1 -> Contactor
- R1 -> Resistencia / *Heater*

**DESA 0041**

## 1.11. Repuestos

## 1.11. Spare parts

YE.DRDES 00110 Resistencia 110 V / *Heater for 110 V*

YE.DRDES 00220 Resistencia 230 V (1 und p/ DESA 0041, 2 unds p/ DESA 0081)  
*Heater for 230 V (1 unit f/ DESA 0041, 2 units f/ DESA 0081)*

YP.DJUSD 00035 Junta silicona para DESA 0041 / *Silicone gasket for DESA 0041*

YP.DJUSD 00075 Junta silicona para DESA 0081 / *Silicone gasket for DESA 0081*

YE.DSONI 00001 Interruptor de nivel para DESA / *Level switch for DESA*

YE.DTERM 00001 Termostato para DESA / *Thermostat for DESA*