

POBEL

MATERIAL GENERAL DE LABORATORIO
LABORATORY GLASSWARE & EQUIPMENT



TERMOMETRO DIGITAL, A BATERIAS DE 1,5 V. (ag13) X 1 - INTERRUPTOR ON/OFF / TECLA °C/°F / TECLA HOLD PARA FIJAR TEMPERATURA / TECLA MAX/MIN.

DIGITAL THERMOMETER - POWER SUPPLY 1.5V(AG13) X 1 - ON/OFF KEY SWITCHING - °C/°F DISPLAY MODE, HOLD KEY AND MAX/MIN KEY

REFERENCIA CAT. NUMBER	RANGO / RANGE	PRECISION / PRECISION
10TEDI0010	-50°C/+300°C / - 58°F + 572°F	± 1°C/ ±1°F



RELOJ AVISADOR MECANICO PROGRAMABLE DE 0-60 MINUTOS. SEÑAL ACÚSTICA DIV. 1 MIN. TIMER 0 TO 60 MINUTES, AUDIBLE ALARM. DIV. 1 MIN

REFERENCIA CAT. NUMBER
10RAEL0010



TERMOMETRO DE LABORATORIO, CON FRANJA AMARILLA, Y ANILLA PARA COLGAR LABORATORY THERMOMETER, WITH YELLOW BRAND AN RING TO HAND

REFERENCIA CAT. NUMBER	GRADUACION GRADUATION	DIVISIONES DIVISION	LIQUIDO /LIQUID
10TERM0110	-35°C/+110°C	1/1°C	ROJO
10TERM0250	-10°C/+250°C	1/1°C	ROJO



PINZA EN ABS FORCEPS IN ABS

REFERENCIA CAT. NUMBER
10PINZ0050



ESPATULA/AGITADOR DE PP. SPATULE/ROD IN PP

REFERENCIA CAT. NUMBER
10ESPA0010




MICRO HOMOGENEIZADOR PP PARA TUBOS DE 0,5 Y 1,5 ML MICRO PESTLE FOR MICROCENTRIFUGE TUBES 0,5 AND 1,5 ML. IN PP

REFERENCIA CAT. NUMBER
10MIMA0010




TUBOS DE CENTRIFUGA FONDO REDONDO CON TAPA A PRESION DE PP. AUTOCLAVABLES POLIPROPYLENE TUBES ROUND BOTTON WITH CAP

REFERENCIA CAT. NUMBER	CAPACIDAD CAPACITY	PAQUETE BAG QUANTITY	MEDIDAS / MEASURES
10TCFR0015	15 ML.	100	17x100 MM.
10TCFR0050	50 ML.	100	28X102 MM.



ESTUCHE PARA PORTA-OBJETOS SLIDE BOX

REFERENCIA CAT. NUMBER	CANTIDAD
10ESPO0050	PARA / FOR 50
10ESPO0150	PARA / FOR 100



MATRACES AFORADOS EN PP. CON TAPON A ROSCA , AUTOCLAVABLES VOLUMETRIC FLASK IN PP WITHN SCREW CAPS. AUTOCLAVABLE.

REFERENCIA CAT. NUMBER	CAPACIDAD CAPACITY
10MAPR0025	25 ML.
10MAPR0050	50 ML.
10MAPR0100	100 ML.
10MAPR0250	250 ML.
10MAPR0500	500 ML.